

User Guidelines for

GEL-MA INX<sup>®</sup> X100



## General Information

### Storage

GEL-MA INX should be stored in a fridge at 4 °C until ready to use. Protect it from light. Expiry date of the product is indicated on the sealed pouch. The product can be stored for a maximum of 3 months after opening and should be consumed before the expiry date.

### Intended Use

Research use only. This product is not intended for use in diagnostic or therapeutic procedures.

### Safety Information

For more information, please refer to the material safety data sheet.

## User Guidelines

### Preparation



GEL-MA INX X100 was produced under sterile conditions. To ensure optimal performance and prevent contamination, it is recommended to handle this product in a **sterile environment**.

1. Warm up the cartridge at room temperature for 20 min.
2. Pre-heat the printhead of the 3D printer at 27 °C.
3. Remove the end-cap and tip-cap of the cartridge and attach a printing nozzle. Insert the cartridge in the pre-heated printhead (27 °C), and warm up for at least 20 min before printing. (See processing guidelines for recommended nozzle types).



For an optimal printing performance, the use of a metal conical **nozzle insulator** is required. Insert the nozzle tip in the insulator as shown in the image.



4. Calibrate the printhead and start printing using the suggested printing parameters (See processing guidelines)

## Processing

### a) Recommended parameters for pneumatic-based printers

Recommended processing parameters for a pneumatic-based 3D printer are listed below.

	<b>22G Nozzle</b>	<b>25G Nozzle</b>	<b>27G Nozzle</b>
<b>Nozzle geometry</b>	Conical	Conical	Conical
<b>Printhead Temperature</b>	27 ± 0.5 °C	27 ± 0.5 °C	27 ± 0.5 °C
<b>Printbed Temperature</b>	16 ± 1 °C	16 ± 1 °C	16 ± 1 °C
<b>Pressure</b>	60 ± 20 kPa	70 ± 20 kPa	80 ± 20 kPa
<b>Infill Speed</b>	5 ± 1 mm/s	5 ± 1 mm/s	5 ± 1 mm/s
<b>Layer Height</b>	0.21 ± 0.02 mm	0.14 ± 0.02 mm	0.11 ± 0.02 mm

### b) Recommended parameters for piston-based printers

Recommended processing parameters for a piston-based 3D printer are listed below.


	<b>22G Nozzle</b>	<b>25G Nozzle</b>	<b>27G Nozzle</b>
<b>Nozzle geometry</b>	Conical	Conical	Conical
<b>Printhead Temperature</b>	27 ± 0.5 °C	27 ± 0.5 °C	27 ± 0.5 °C
<b>Printbed Temperature</b>	16 ± 1 °C	16 ± 1 °C	16 ± 1 °C
<b>Flow speed</b>	3 ± 1 mm/s	3 ± 1 mm/s	3 ± 1 mm/s
<b>Infill Speed</b>	5 ± 1 mm/s	5 ± 1 mm/s	5 ± 1 mm/s
<b>Layer Height</b>	0.21 ± 0.02 mm	0.14 ± 0.02 mm	0.11 ± 0.02 mm

 The printing parameters have been validated for printing a cube with dimensions 15 x 15 mm (W x L) using 3 ml cartridges.

 Use of rectilinear infill pattern is recommended.

**! Photo-crosslinking:** During printing, structure should be irradiated with light ( $\lambda$ : 365 nm or 405 nm, Dose: 70 mJ/cm<sup>2</sup>) after every five layers. This step is required for partial crosslinking of the structure for a better shape retention. After completion of printing, the final structure should be promptly placed under UV light for complete crosslinking. (Recommended parameters for post-printing photo-crosslinking:  $\lambda$ : 365 or 405 nm, Dose: 10000 mJ/cm<sup>2</sup>)

For photo-crosslinking kinetics of GEL-MA INX at two different wavelengths, see Figure 1.



**Figure 1. Storage ( $G'$ , solid lines) and loss ( $G''$ , dashed lines) moduli of GEL-MA INX X100 as a function of irradiation dose at 365 nm and 400-500 nm wavelengths**

## Cell Culture

### a) Cell Seeding

The scaffolds can be readily seeded with cells after overnight incubation in cell culture media.

### b) Cell Encapsulation



Preheat GEL-MA INX cartridge at 37 °C.



Connect GEL-MA INX cartridge to a syringe via a female luer-to-luer adapter.



Inject the desired amount of GEL-MA INX into the syringe.



Connect the second syringe containing the cell suspension via a luer-to-luer adapter. Shift the hydrogel back and forth between the syringes until the cells are homogeneously mixed in.



Inject cell-laden GEL-MA INX back into the cartridge via a luer-to-luer adapter.