

Specification Sheet

Product Name	BIORES INX X100
Product Code	DLP-101
Appearance	Yellow translucent liquid (20 °C)
Total Degree of Functionalization	90 – 100 %
Viscosity	5 – 20 mPa.s
Tested using rheometer at 20 °C	
Storage Modulus after UV curing	20 – 40 kPa
Tested using rheometer at 20 °C, (wavelength: 400 – 480 nm, Intensity: 20 mW cm ⁻²)	

Related Documents

Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

Intended Use

For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.