

Specification Sheet

Product Name	BIORES INX X100
Product Code	DLP-101
Appearance	Yellow translucent liquid (20 °C)
Total Degree of Functionalization	90 – 100 %
Viscosity	10 – 30 mPa.s
Tested using rheometer at 20 °C	
Storage Modulus after UV curing	25 – 50 kPa
Tested using rheometer at 20 °C, (wavelength: 400 – 480 nm, Intensity: 20 mW cm ⁻²)	

Related Documents	Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.
--------------------------	---

Intended Use	For research use only. Not intended for human use <i>in vivo</i> . Not to be used in diagnostic, therapeutic or clinical applications.
---------------------	--