

Specification Sheet

Product Name DEGRES INX X100

Product Code DLP-301

Appearance Yellow liquid

Viscosity 0.05 – 0.60 Pa.s

Tested using rheometer at 20 °C

Storage Modulus after UV curing 50 – 500 kPa

Tested using rheometer at 20 °C,

(wavelength: 400 – 480 nm, Intensity: 20 mW cm⁻²)

Related Documents

Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

Intended Use

For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.