

Specification Sheet DEGRAD INX

Version 1.0 September 2021

Specification Sheet

Product Name DEGRAD INX

Product Code MPL-301

Appearance Yellow-orange liquid

Viscosity 0.5 - 3 Pa.s

Tested using Anton Paar Physica MCR-301 Rheometer at 20°C

Related Documents Safety Data Sheet and Certificate of

Analysis are provided by the manufacturer.

Additional information concerning material

handling can be found in the Application

Protocol.

Intended Use For research use only. Not intended for

human use in vivo. Not to be used in

diagnostic, therapeutic or clinical

applications.

Confidential 1