

## Specification Sheet

Product Name	HYDROTECH INX X100
Product Code	MPL-201
Appearance	Yellow - orange liquid
Viscosity	0.5 – 5 Pa.s

Tested using Anton Paar Physica MCR-301  
Rheometer at 20 °C

### Related Documents

Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

### Intended Use

For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.