

## Specification Sheet

<b>Product Name</b>	HYDROTECH INX X200
<b>Product Code</b>	MPL-202
<b>Appearance</b>	Orange liquid at 20 °C
<b>Viscosity</b>	1 – 10 Pa.s
<small>Tested using rheometer at 20 °C</small>	
<b>Storage Modulus</b>	500 – 1200 kPa
<small>Tested using rheometer at 20 °C, 300-390 nm, 20 mW cm<sup>-2</sup></small>	
<b>Related Documents</b>	Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.
<b>Intended Use</b>	For research use only. Not intended for human use <i>in vivo</i> . Not to be used in diagnostic, therapeutic or clinical applications.