

Specification Sheet

Product Name	HYDROTECH INX X100
Product Code	MPL-201
Appearance	Yellow - orange liquid
Viscosity	0.5 – 5 Pa.s

Tested using Anton Paar Physica MCR-301 Rheometer at
20 °C

Related Documents

Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

Intended Use

For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.