



X-PURE® GeIDAT 120P30

Technical grade

Product description

X-Pure® GeIDAT is Gelatin Desaminotyrosin. The precursor gelatin is an acid extracted pig skin gelatin. X-pure® GeIDAT is specifically designed for applications requiring intensely purified cross-linkable biomaterials.

X-Pure Technical Grade products are physically equivalent (DoM & MW) to the research and GMP grades. Rousselot provides no warranties with regards to the material purity or other non-physical characteristics for the Technical grade material.

Technical grade materials are intended solely for physical material characterization (e.g. optimization of hydrogel strength).

X-Pure® GeIDAT appears as white to off-white lyophilized powder.

Physical/chemical limits

Standard parameters	Typical values	Test Method
Degree of modification	20 – 40%	Rousselot
Molecular weight	100 – 140 kDa	Rousselot

Packaging

PET or PETG flask in an aluminum bag. Quantity 50g and 10g.

Storage information

Store away from heat and moisture.

Caution: Close the recipient immediately after use.

Rousselot Biomedical

As the most recent strategic segment within Rousselot, we have drawn upon Rousselot's 130+ years of worldwide expertise and proven track record of pharmaceutical gelatins and collagens to develop innovative ranges of purified, modified and non-modified gelatins and collagens for biomedical applications. Offering unique advantages to assure performance, quality and safety from bench to clinic, Rousselot® X-Pure® and Rousselot® Quali-Pure® are backed by strong scientific data and on-going research. Rousselot Biomedical is committed to support end-to-end partnerships to help "Improving patients' lives".

Rousselot is Darling Ingredients' Health Brand.



Supplier

Rousselot BV
Meulestedekaai 81 – 9000 Gent – Belgium
Tel: + 32 (0)9 255 18 18
Biomedical@rousselot.com

*Non contractual document, specifications presented can be changed without prior notice. Ensure you use latest version /
Reference PDS-BIOMED-10011 Revision: 2 Date: Feb 19, 2024 © Rousselot BV - Pictures: © Rousselot /*