

Specification Sheet

Product Name READYPCL INX X100

Product Code VOL-301

Appearance Clear liquid

After heating at 50 °C

Viscosity 5 – 20 Pa.s

Tested using rheometer at 20 °C

Storage Modulus after UV curing 1 – 5 MPa

Tested using rheometer at 20 °C,
(wavelength: 400 – 500 nm, Intensity: 20 mW cm⁻²)

Refractive Index 1.46 – 1.48

Tested at 589 nm at 20 °C

Related Documents Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

Intended Use For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.