



## Specification Sheet

<b>Product Name</b>	GEL-MA INX X100
<b>Product Code</b>	DEP-101
<b>Appearance</b>	Transparent gel
<b>pH</b>	6.5 – 8.5
<b>Total degree of functionalization</b>	70 - 80%

<sup>1</sup>H NMR performed in D<sub>2</sub>O at 40°C, acquired with a spectral width of 8013 Hz, or 16 ppm, averaged over 64 scans using 64K time domain points.

### Related Documents

Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

### Intended Use

For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.