



Specification Sheet

Product Name	STABLE INX X100
Product Code	DEP-401
Appearance	Transparent gel
pH	7.0 – 9.0
Viscosity at high shear	0.5 – 5 Pa·s

Tested using Anton Paar Physica MCR-301 Rheometer @ 25°C

Viscosity at low shear	1000 – 7000 Pa·s
-------------------------------	------------------

Tested using Anton Paar Physica MCR-301 Rheometer @ 25°C

Storage modulus after UV crosslinking	0.25 – 1.00 MPa
--	-----------------

Tested using Anton Paar Physica MCR-301 Rheometer @ 25°C.

Related Documents

Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

Intended Use

For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.