

Specification Sheet

Product Name SOLID INX X100

Product Code DEP-501

Appearance Off-white solid

Melting Temperature 40 – 60°C

DSC analysis with the use of TA instruments, (Q2000 DSC).

Crystallinity 30 – 40 %

DSC analysis with the use of TA instruments, (Q2000 DSC).

Viscosity 10 - 30 Pa·s

Tested using Anton Paar Physica MCR-301 Rheometer
@ 60°C.

Related Documents

Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

Intended Use

For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.