

## Specification Sheet

**Product Name** X SUPPORT INX

**Product Code** DEP-601

**Appearance** Transparent gel

**pH** 9 – 11

**Viscosity at high shear** 0.1 - 1 Pa·s

Tested using Anton Paar Physica MCR-301 Rheometer @  
25°C

**Viscosity at low shear** 3000 – 7000 Pa·s

Tested using Anton Paar Physica MCR-301 Rheometer @  
25°C

### Related Documents

Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

### Intended Use

For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.