

Specification Sheet

Product Name	X-PURE GEL-MA INX X100
Product Code	DEP-111
Appearance	Transparent gel
pH	6.5 – 8.5
Viscosity	0.01 – 500 Pa.s
Tested with rheometer at 27 °C	
Storage Modulus	5 – 15 kPa
Tested with rheometer at 27 °C during UV irradiation	

Related Documents

Safety Data Sheet and Certificate of Analysis are provided by the manufacturer. Additional information concerning material handling can be found in the Application Protocol.

Intended Use

For research use only. Not intended for human use *in vivo*. Not to be used in diagnostic, therapeutic or clinical applications.